

# MONTHLY DISPOSABLE CONTACT LENSES

## FREEFORM STABILIZATION

Avanti's freeform stabilization utilizes Peripheral Balancing Zones (PBZ) that interact with the top and bottom lids. This interaction creates stabilizing forces which allow for rapid stabilization, usually within 20 seconds.

## MULTIFOCAL OPTIONS

Avanti® is available in Centre Distance and Centre Near designs with the additional option to optimize for the dominant eye.

Tell us your patient's needs and our customer service team will advise the most suitable combination of options.

## ADVANTAGES

- Wide parameter range
- Stable vision through gaze excursions
- Does not mislocate during wear
- Stability maintained when head is tilted
- Reduction in overall thickness
- Increased comfort

THE CONVENIENCE OF A MONTHLY WITH THE FLEXIBILITY OF CUSTOM MADE



MAKE SURE YOUR PATIENTS NEVER MISS THE ACTION

CALL NOW FOR YOUR FREE TRIALS

 **0800 585115**

Tel +44 (0) 1252 381112 • [info@ultravision.co.uk](mailto:info@ultravision.co.uk) • [www.ultravision.co.uk](http://www.ultravision.co.uk)

  
**avanti**<sup>®</sup>  
MONTHLY DISPOSABLE SOFT CONTACT LENSES

Avanti® is a registered trade mark of UltraVision International Limited.

  
INNOVATION IN PRACTICE

  
**avanti**<sup>®</sup>  
MONTHLY DISPOSABLE SOFT CONTACT LENSES

FAST  
FOCUS  
FINE  
TUNED



# MAKE SURE YOUR PATIENTS NEVER MISS THE ACTION

Avanti's unique design allows for rapid stabilization ensuring crisp, clear vision at all times. Gaze excursions and tilting your head is no longer a problem as lenses do not mislocate.

**A truly flexible Multifocal system** combines Centre Near and Centre Distance designs with additional options to optimize vision for the dominant eye.

The superb stability of Avanti® and its comprehensive range of visual options ensure your patients are always ahead of the action.

  
**avanti**®  
MONTHLY DISPOSABLE SOFT CONTACT LENSES



## FAST FOCUS FINE TUNED

### PARAMETER RANGE

Base Curve	8.00 8.30 8.60 8.90 9.20 mm
Diameter	14.00 14.50 15.00 mm
Sphere	+/- 20.00 DS
Cyl	-0.75 to -8.00 DC
Axis	1°- 180° (1° steps)
Multifocal	Centre Distance & Centre Near
Add	+0.50 to +4.00 DS

### MATERIAL

Definitive 74% SiH (Blue handling tint)  
Dk 60  
UV Protection (Class 2)

### EASY TO ORDER

Call or e-mail with HVID, K readings and Spectacle Rx to receive your diagnostic lenses.